

## Which Z2 Is Right For You?

### Z2 CPAP

The Z2® CPAP machine is an ultra-small, light weight, full featured fixed-pressure cpap machine. It is the ultimate fixed pressure cpap for travelers and people with an active lifestyle.


### Z2 Auto

The Z2 Auto is the same size and weight as the fixed pressure CPAP making it perfect for travel. The therapy pressure adjusts on a breath by breath basis which results in a lower overall mean pressure than the fixed pressure CPAP.

	Z2 CPAP	Z2 Auto
Fixed Pressure Mode	✗	✗
Auto Adjusting Pressure Mode		✗
Auto Stop Feature	✗	
Auto Start / Stop Feature		✗
Approved for In-Flight Use	✗	✗
PC Data Capable	✗	✗
Mobile Data Capable		✗
Pressure Range 4-20 cm H <sub>2</sub> O	✗	✗
Ramp 0-45 minutes	✗	✗
Z-Breathe® Comfort Feature	✗	✗
Auto Adjusting Voltage for World Travel (100-240v)	✗	✗
Qtube Included	✗	✗
No Special Mask Required	✗	✗

## Check Your Treatment Information

The backlit LCD Screen shows your treatment settings on the device. You can easily adjust your comfort settings directly on your device.



1. Treatment pressure display
2. Ramp progress indicator (if ramp set)
3. Ramp time remaining (if ramp set)
4. Battery status (if present)
5. AC Power indicator (if connected)
6. General status display (date, mode, Auto Start/Stop)
7. Fault indicator (if fault occurs)

If you travel for work or just for fun  
ask about the Z2 today.



[CLICK TO LEARN MORE](#)

**AUTHORISED DISTRIBUTOR**

**HITWIN HEALTHCARE PVT LTD**

7/10, Jaganathan Nagar 2nd Main Road, Arumbakkam,  
Chennai - 600106. India.

**Hitwin**  
HEALTHCARE

**BREAS**  
A GE Healthcare Company

Breas Medical/2019/BOS-5AL-NAR-000337-v1.0

**Z2 auto™**  
by Breas  
Ultra-Small, Travel CPAP  
**New!**



Quieter Motor  
Auto Start / Stop  
No Special Mask Required

## Take Your Adventures Off the Grid

Whether you need the fixed pressure Z2<sup>®</sup> CPAP or the auto-adjusting Z2 Auto,<sup>®</sup> you can take your CPAP anywhere you want to go. Just add the optional PowerShell<sup>®</sup> with Extended Life Battery to provide battery power for a full night's sleep on a single charge.

It's great for camping or any where access to power is limited.



Shown Actual Size  
6.48" x 3.30" x 2.02"




## Travel Anywhere in the World

Both Z2s are FAA Compliant for in-flight use. They come with auto-adjusting voltage power supplies that let you plug into any outlet in the world with just an added plug adapter.

## Check Your Sleep Data

Download your Z2 Auto<sup>®</sup> sleep data directly to your iOS or Android phone or tablet with the free Nitelog<sup>®</sup> mobile app.

Just download the app, connect to the Z2 via Bluetooth and you're ready to go.



**Hitwin**  
HEALTHCARE

Respiratory • Sleep Therapy • Cardiology