

INTEGRATED HEATED TUBE



LH1

Integrated heated tube



Integrated Heating Design

Integrated heating solution minimizes the condensate and improves the compliance.

Compatible

It is designed with 22mm cuffs so that it is compatible with most of the CPAP devices.

Specifications

Tube	Length	Power supply	Compliance standards	Resistance to flow at rated flow
LH1	1.8m	24 V / 18 W	60 hPa: < 5 ml/hPa	30 l/min: < 0.06 hPa/l/min



Kink Resistant

Scientific spiral interval design makes the tubing more flexible and kink resistant.

