

Wardcare Solutions

Fowler (Back and Knee / Leg adjustable) bed with conventional head and foot boards and side railing



Fowler (Back and Knee / Leg adjustable) bed with stainless steel laminated head and foot boards and side railing

Fowler (Back and Knee / Leg adjustable) bed with polymer moulded head and foot boards, and side railing



BX4000 Plus Fowler Bed

The BX4000 plus fowler bed has the facility to raise the back and upper leg sections, which are operated by the crank mechanism provided at the foot end of the bed. The raising section allows the nursing staff to change positions for various activities and clinical positions with relative ease.

*The product images shown are for illustrative purpose only and not be an exact representation of the product. Accessories shown in the product image may not be part of the standard configuration.

Technical Specifications

Bed Frame Dimension*	2080 mm x 920 mm, 82" x 36" (L x W)
Mattress Platform*	1970 mm x 820 mm, 78" x 32" (L x W)
Buffer to Buffer * (when castors provided)	2130 mm x 970 mm, 84" x 38" (L x W)
Backrest Angle**	0° - 65°
Knee Rest Angle**	0° - 50°
Safe Working Load	135 kgs
Patient Weight Bearing Capacity^	100 kgs

Standard Features

- Raised backrest
- Raised kneerest
- DVT end platform/calf elevation
- Four section metallic perforated top
- Pre-treated and epoxy powder coated MS parts
- Provision for IV Pole
- Removable handle
- Easily removable head and foot boards
- Supplied in SKD condition

Safety Standard

General Safety for Medical Bed IEC 60601-2-52

Add-On Accessories

Railings and Legs

- Collapsible type epoxy powder coated vertical safety side railing (pair)
- H' type legs 910 mm (36") width & 600 mm (24") height with rubber shoes
- H' type legs 910 mm (36") width & 450 mm (18") height with rubber shoes
- H' type legs 910 mm (36") width & 600 mm (24") height with 125 mm high grade synthetic castors, two with brake and two without brake
- H' type legs 910 mm (36") width & 450 mm (18") height with 125 mm high grade synthetic castors, two with brake and two without brake

- Four section mattress with 4" thick PU foam of 40 density covered with PVC
- Telescopic IV pole with 2 hooks
- Telescopic IV pole with 4 hooks
- Orthopaedic traction pulley for trauma cases
- Fixed polymer moulded handle (set of 2)
- Removable handle
- Moulded chart holder
- Mattress arrestor (set of 2)
- 4" safety rubber buffer at four corners
- Urine bag holder

Head and Foot Boards

- Polymer moulded head and foot boards
- SS laminated head and foot boards
- Conventional head and foot boards

*± 10 mm Engineering Variation and ** ± 5° Engineering Variation ^ When Bed is configured with standard accessories
Note: Option of fixed polymer moulded handle available in lieu of removable handle.

