

Hospital Furniture

KRAFT™
Crafted for Care & Comfort




Hitwin
HEALTHCARE

Hitwin Healthcare Pvt. Ltd.

7/10, Jaganathan Nagar 2nd Main Road,
Arumbakkam, Chennai, India - 600106.

Ph: 91 44 4684 5555 (20 Lines), 91 44 4233 3313

Mob: 91 93441 66778




KRAFT 100
MOTORIZED ICU BED
5 FUNCTION (EXCEL PLUS)




KRAFT 101
MOTORIZED ICU BED
5 FUNCTION (EXCEL)

TECHNICAL SPECIFICATION FOR 5 FUNCTION ICU BEDS

- Size(L×W×H): 2150 mm × 965 mm × 485 to 710 mm (height adjustable)
- The motorized bed comes with 4 actuators with control box & handset
- Hydraulically pressed four section bed boards in 18 SWG CRCA sheet
- Head & Leg Bows set made of high quality engineering plastic (set of 2)
- Tuck away type ABS polymer side railings (set of 4)
- Provision for IV Pole, fitted with drainage bag holder
- Mobile on 125 mm castors 4 nos with central locking system
- Epoxy powder coating with 8 tank process for scratch resistance
- Optional (at extra cost):
 - Nursing panel
 - Battery backup for Electric Bed



KRAFT 102
MOTORIZED ICU BED
5 FUNCTION (PREMIUM)



KRAFT 103
MANUAL ICU BED
5 FUNCTION (DELUXE)

TECHNICAL SPECIFICATION FOR 5 FUNCTION ICU BEDS

- Size(L×W×H): 2150 mm × 965 mm × 485 to 710 mm (height adjustable)
- The motorized bed comes with 4 actuators with control box & handset. Manual bed comes with 4 cranks
- Hydraulically pressed four section bed boards in 18 SWG CRCA sheet
- Head & Leg Bows set made of high quality engineering plastic (set of 2)
- Aluminium collapsible side railings with locking system (set of 2)
- Provision for IV Pole, fitted with drainage bag holder
- Mobile on 125 mm castors 4 nos with all brakes
- Epoxy powder coating with 8 tank process for scratch resistance
- Optional (at extra cost):
 - Nursing panel
 - Battery backup for Electric Bed
 - Castors with central locking mechanism



KRAFT 121
MOTORIZED ICU BED
3 FUNCTION (PREMIUM)



KRAFT 122
MANUAL ICU BED
3 FUNCTION (DELUXE)

TECHNICAL SPECIFICATION FOR 3 FUNCTION ICU BEDS

- Size(L×W×H): 2150 mm × 965 mm × 485 to 710 mm (height adjustable)
- The motorized bed comes with 3 actuators with control box & handset. Manual bed comes with 3 cranks
- Hydraulically pressed four section bed boards in 18 SWG CRCA sheet
- Head & Leg Bows set made of high quality engineering plastic (set of 2)
- Aluminium collapsible side railings with locking system (set of 2)
- Provision for IV Pole, fitted with drainage bag holder
- Mobile on 125 mm castors 4 nos with all brakes
- Epoxy powder coating with 8 tank process for scratch resistance
- Optional (at extra cost):
 - Nursing panel
 - Battery backup for Electric Bed
 - Castors with central locking mechanism





KRAFT 131
MOTORIZED FOWLER BED
2 FUNCTION (PREMIUM)



KRAFT 132
MANUAL FOWLER BED
2 FUNCTION (DELUXE)



KRAFT 133
MANUAL FOWLER BED
2 FUNCTION (ECONOMY)

TECHNICAL SPECIFICATION FOR FOWLER (2 FUNCTION) BEDS

- Size (L×W×H): 2050-2100 mm × 900-960 mm × 500 mm
- The motorized bed comes with 2 actuators with control box & handset. Manual bed with 2 cranks
- Four section bed boards in 18 SWG CRCA sheet
- Head & Leg Bows set made of high quality engineering plastic (set of 2)
- Provision for IV Pole, fitted with drainage bag holder
- Epoxy powder coating with 8 tank process for scratch resistance
- Optional :
 - Aluminium collapsible side railings with locking system (set of 2)
 - Castors 125mm dia - 4 nos all with brakes



KRAFT 141
MOTORIZED SEMI FOWLER
BACKREST BED
1 FUNCTION (PREMIUM)



KRAFT 142
MANUAL SEMI FOWLER
BACKREST BED
1 FUNCTION (DELUXE)

TECHNICAL SPECIFICATION FOR SEMI - FOWLER (1 FUNCTION) BEDS

- Size (L×W×H): 2100 mm × 900 mm × 500 mm
- The motorized bed comes with 1 actuator with control box & handset. Manual bed comes with 1 crank
- Two section bed boards in 18 SWG CRCA sheet
- Head & Leg Bows set made of high quality engineering plastic (set of 2)
- Provision for IV Pole, fitted with drainage bag holder
- Epoxy powder coating with 8 tank process for scratch resistance
- Optional :
 - Aluminium collapsible side railings with locking system (set of 2)
 - Castors 125mm dia - 4 nos all with brakes



KRAFT 143
MANUAL SEMI FOWLER
BACKREST BED
1 FUNCTION (ECONOMY)



KRAFT 144
MANUAL SEMI FOWLER
BACKREST BED
1 FUNCTION (CLASSIC)

TECHNICAL SPECIFICATION FOR SEMI - FOWLER (1 FUNCTION) BEDS

- Size (L×W×H): 2000 mm × 900 mm × 500 mm
- Two section bed boards in 20 SWG CRCA sheet
- The backrest function is operated with crank
- Head & Leg Bows set made of high quality engineering plastic (set of 2)
- Provision for IV Pole, fitted with drainage bag holder
- Epoxy powder coating with 8 tank process for scratch resistance
- Optional :
 - Aluminium collapsible side railings with locking system (set of 2)
 - M.S side rails (set of 2) for KRAFT 144
 - Castors 125 mm dia - 4 nos all with brakes



Ward / OT Furniture



KRAFT 151
PLAIN BED CLASSIC
SIZE : L 1900mm X W 900 mm X H 560mm



KRAFT 152
ATTENDANT BED
SIZE : L 1830mm X W 610mm X H 380mm



KRAFT 153
PEDIATRIC BED CLASSIC
SIZE : L 1370mm X W 625mm X H 560mm



KRAFT 154
ATTENDANT BED CUM CHAIR



KRAFT 201
EXAMINATION COUCH
SIZE : L 1830mm X W 575 mm X H 850mm



KRAFT 202
EXAMINATION COUCH WITH BACKREST
SIZE : L 1830mm X W 575 mm X H 850mm




KRAFT 203
EXAMINATION TABLE S.S.
SIZE : L 1830mm X W 537mm X H 810mm



KRAFT 204
SIMPLE EXAMINATION COUCH
SIZE : L 1830mm X W 575mm X H 850mm



KRAFT 205
BLOOD DONATION COUCH



KRAFT 301
MOTORIZED OBSTETRIC LABOR & RECOVERY BED - 3 FUNCTION



KRAFT 302
TELESCOPIC OBSTETRIC & GYNAEC TABLE - S.S.
SIZE : L 1830mm X W 760 mm X H 760mm



KRAFT 303
SIMPLE GYNAEC DELIVERY TABLE S.S./M.S.
SIZE : L 1830mm X W 760mm X H 760mm



KRAFT 304
MOTORIZED GYNAEC OT TABLE



KRAFT 400
TRAUMA CARE RECOVERY TROLLEY
By Screw Mecaahnism



KRAFT 401
HYDRAULIC TRAUMA CARE
RECOVERY TROLLEY



KRAFT 402
STRETCHER ON TROLLEY - SS
SIZE : L 2030mm X W 610mm X H 870mm



KRAFT 403
ABS PATIENT TRANSFER TROLLEY



KRAFT 404
ALUMINIUM FOLDING STRETCHER



KRAFT 500 A
ABS EMERGENCY TROLLEY



KRAFT 500 B
ABS ANESTHESIA TROLLEY



KRAFT 501
EMERGENCY CRASH CART - S.S



KRAFT 502
BABY TROLLEY - S.S
SIZE : L 750mm X W 455mm X H 1000mm



KRAFT 503
MONITOR TROLLEY M.S



KRAFT 504
WASTE CARRYING TROLLEY S.S



KRAFT 505
BOWL STAND SINGLE - S.S
SIZE : H 840 mm X D 355 mm



KRAFT 506
BOWL STAND DOUBLE - S.S
SIZE : H 840 mm X D 355 mm



KRAFT 507
REVOLVING STOOL - S.S
SIZE : W 290mm X H 425mm to 675mm



KRAFT 508
SOILED LINEN TROLLEY BIG - S.S
SIZE : L 790mm W 560mm X H 1030mm



KRAFT 509
BED SIDE STOOL WITH S.S TOP
SIZE : L 405mm X W 405 mm X H 455mm



KRAFT 510 / 511
SINGLE STEP STOOL - S.S / M.S
SIZE : L 450mm X B 225mm X H 225mm
DOUBLE STEP STOOL - S.S / M.S
SIZE : L 450mm X B 400mm X H 400mm



KRAFT 512
SALINE STAND - S.S



KRAFT 601
DRESSING TROLLEY
Available in various sizes



KRAFT 602
INSTRUMENT TROLLEY - S.S
Available in various sizes



KRAFT 603
BED SIDE TROLLEY - S.S
Available in various sizes



KRAFT 604
INSTRUMENT TROLLEY - S.S
Available in various sizes

S.S. TROLLEY - VARIOUS SIZES

L X W X H
16 X 16 X 35
18 X 18 X 35
24 X 18 X 35
27 X 18 X 35
30 X 20 X 35
36 X 20 X 35
36 X 24 X 35
42 X 24 X 35



KRAFT 605
OVER BED TABLE - S.S
SIZE : L 1115mm X W 445mm X H 965mm



KRAFT 606
MAYO'S TROLLEY - S.S
SIZE : L 560mm X W 445mm X H 810 to 1270mm



KRAFT 701
BED SIDE LOCKER LAMINATED
SIZE : L 380mm X W 380mm X H 820mm



KRAFT 702
BED SIDE LOCKER ABS
SIZE : L 475mm X W 475mm X H 820mm



KRAFT 703 A
ADJUSTABLE BED SIDE TABLE -
GAS SPRING



KRAFT 703 B
ADJUSTABLE BED SIDE TABLE
WITH KNOR



KRAFT 704
NON FOLDING WHEEL CHAIR S.S



KRAFT 704 A
FOLDING WHEEL CHAIR



KRAFT 705
BED SIDE SCREEN S.S / M.S
SIZE : H 1675mm X W 2440mm



KRAFT 706
CYLINDER TROLLEY
S.S / M.S 10 LTRS



KRAFT 801
PP HEAD & LEG BOWS (ECONOMY)
approx wt: 5.500 kgs per set



KRAFT 802
PP HEAD & LEG BOWS (DELUXE)
approx wt: 7.500 kgs per set



KRAFT 803
PP HEAD & LEG BOWS (PREMIUM)



KRAFT 804
HEAD & LEG BOWS (EXCEL)
approx wt: 8.400 kgs per set




KRAFT 805
HEAD & LEG BOWS (EXCEL PLUS)
approx wt: 9.800 kgs per set




KRAFT 806
TUCK AWAY PP SIDE RAILS
(SET OF 4)



KRAFT 807
TUCK AWAY PP SIDE RAILS
(SET OF 4)



KRAFT 808
TUCK AWAY PP SIDE RAILS
(SET OF 2)



KRAFT 809
STEEL COLLAPSIBLE SIDE RAILS
(SET OF 2)
approx wt: 5.500 kgs per set



KRAFT 810
ALUMINIUM HEAVY DUTY COLLAPSIBLE
SIDE RAILS (SET OF 2)
approx wt: 6.100 kgs per set




KRAFT 811
SINGLE FOLD FOAM MATTRESS



KRAFT 812
MULTI FOLD FOAM MATTRESS



KRAFT 813
CRASH CART ACCESSORIES
2 CABINET & 6 UTILITY BINS



KRAFT 814
BED CRANK - 48"



KRAFT 815
ABS EXPANDABLE FOOD TRAY



KRAFT 816
CENTRAL LOCKING CASTORS
125MM




KRAFT 817
TWIN CASTORS
125MM



KRAFT 818
HEAVY CASTORS
125MM



KRAFT 819
ECONOMY CASTORS
125MM



KRAFT 820
ABS CHART HOLDER




KRAFT 821
BED SIDE TABLE TOP &
GAS SPRING



KRAFT 822
I.V. ATTACHMENT
FOR BED

KRAFT™
Crafted for Care & Comfort



Hitwin
HEALTHCARE

Hitwin Healthcare Pvt. Ltd.
7/10, Jaganathan Nagar 2nd Main Road,
Arumbakkam, Chennai, India - 600106.
Ph: 91 44 4684 5555 (20 Lines), 91 44 4233 3313
Mob: 91 93441 66778