




A comfortable mask improves adherence to treatment, leading to faster well-being of the patients



AirFit™ N20



AirFit™ F20



Respiratory • Sleep Therapy • Cardiology

Introducing
AirFit™ 20 series

Perfect Fit. Maximum Comfort.

AirFit™ N20



AirFit™ F20



- ✓ Adapts to your face for perfect fit
- ✓ Maximum seal even in high pressure & movements

TO PLACE ORDERS, CONTACT: Hitwin Healthcare Pvt. Ltd.
 7/10, Jaganathan Nagar 2nd Main Road,
 Arumbakkam, Chennai, India - 600106.
Ph: 91 44 4684 5555 (20 Lines), 91 44 4233 3313
Mob: 91 93441 66778

www.hitwinhealthcare.com



A revolutionary technology

that allows the mask to adapt to your face and provide maximum seal



Respiratory • Sleep Therapy • Cardiology

AirFit™ N20 & AirFit™ F20

fits **>95%** of all patients*

Infinity Seal™ Cushion Technology
Designed to provide consistent support to ensure a confident seal regardless of facial type, therapy pressure or mask movement.

Soft nasal bridge
Shapes to a patient's nasal bridge to foster greater comfort and minimise the likelihood of red marks


Pruned silicone surface
Provides a softer cushion texture around the sensitive nasal region for a comfortable seal and easy fit

Adaptive wings
Accommodate a wide range of facial types by curving to the profile of each wearer

Supportive rim
Designed to deliver stability and security to the adaptive areas of the cushion for a robust seal

Textured chin region
Aims to add increased compliance to the chin area for greater comfort and seal performance during movement

Polished silicone finish
Crafted to seal around the mouth and jawline for a smooth, reliable fit while preventing cushion ride up



Supporting features for all-round comfort & superior experience

Visual freedom

Compact under eye frame, provides full visual freedom

Comfortable

Soft material that's durable and easy to use Velcro tabs

Soft frame

Combination of plastic frame + foam + fabric, provides built-in, no fuss comfort

High strength


Tested upto 400N (40kg) of force applied to frame arm attachments, delivers strength and flexibility.

Upgraded magnetic clips

Colour to colour usability guide with stronger grip

360° movement freedom

360° rotation at elbow-to-frame & elbow-to-hose connection



*Source: RealMed AirFit N20 Internal global fitting study of 159 existing RealMed patients, conducted 12/11/2015 : ID A3637023. RealMed AirFit F20 Internal International fitting study of 80 RealMed and non-RealMed patients, conducted between 08/03/2016 - 22/03/2016. Preliminary patient study - data on file: ID A3637701.