



Nasal Cannula (NC11L / NC11M / NC11S)

The cannula, which is for pediatric use, offers different sizes designed for various-weight patients.

family cares



Nasal Cannula (NC11L / NC11M / NC11S)




Adjustable Slider

Thanks to the slider, secure and comfortable fitting can be achieved without need of an additional strap.

Skin-friendly Stickers

The bilateral stickers are designed under the conception of stability, good air permeability and skin-friendly.



Soft-touching Tubing

Gentle touching on the ears and face so that the friction and pressure on delicate tissues are minimized.

Multiple Patient Groups

Color-coded cannula sizes fulfill requirements from various patient groups.

Part Number	Size	Color	Patient's Weight Range	Flow Range
10410083	Small		2 ~ 8 Kg	2 ~ 10 L/min
10410084	Medium		3 ~ 15 Kg	2 ~ 20 L/min
10410085	Large		12 ~ 22 Kg	2 ~ 25 L/min

