



Nasal Cannula


(NC12L / NC12M / NC12S)

HFNC therapy is able to reduce the dead space effectively while make the moisture, heating and oxygen concentration adjustable so that the comfort and compliance of the patient are improved significantly.

family cares



Nasal Cannula (NC12L / NC12M / NC12S)



Part Number	Size	Color	Flow Range
10410089	Small		10 ~ 70 L/min
10410090	Medium		10 ~ 70 L/min
10410091	Large		10 ~ 85 L/min

Our Leading Distributor:

Hitwin Healthcare Pvt. Ltd.

7/10, Jaganathan Nagar 2nd Main Road,
Arumbakkam, Chennai, India - 600106.

Ph: 91 44 4684 5555 (20 Lines), 91 44 4233 3313

Mob: 91 93441 66778