



**360 Inhaler
Rotation**



**Effective Medicine
Delivery**



**Controlled
Air-Flow**



Easy to Operate



Low Noise



**Long Service
Life**



SAMSON NEBULIZER FOR CHILD & ADULT :

Samson Nebulizer was developed in conjunction with respiratory therapists for the successful treatment of asthma, chronic bronchitis, allergies, and other respiratory disorders

Efficient Nebulization Rate of 0.3ml/min; Small Particle Size of MMAD 3.9µm to reach lower airway for effective nebulization; 12ml Medication Capacity

User-friendly – for a fuss-free application, the nebulizer administers medication at the touch of a button, and a complete set of attachments available

Effective nebulization – breaks down the medication to a median particle, the particle size of 3.9 µm, ideal to reach lower airways, assuring fast and effective treatment

Long-lasting – uncompromising quality that meets the global standards and developed to ensure durability that will last through periods of use

For child and adult with 360-degree nebulization and medication capacity of 12 ml for best respiratory care

The high nebulization rate ensures that it fastens the treatment time and shortens the inhalation time. The MMD is approximately 5 micrometres and the Nebulizer rate is 0.4 millilitre per minute.

Breakage of product, water damage, and wear and tear excluded from the warranty. The warranty covers manufacturing defects in products. The customer needs to send the product to our service center for repair under warranty.

Key Specifications

Usage/Application	Nebulization
Operation Mode	Portable
Size	Compact
Power Consumption	Below 75 W
Noise Level	Less than 60 dBA
Litre Flow Range	6 lpm
Age Group	Adult, Children
Operating Temperature Range	10 to 45 Degree C
Frequency	50 Hz
Voltage	AC 230V

