

# OPTIMO COMPRESSOR NEBULIZER

MODEL: NEB 101

**Hitwin**  
HEALTHCARE

Respiratory • Sleep Therapy • Cardiology



## Highlights

- ~ Compressor Nebulizer
- ~ Durable
- ~ Reliable
- ~ Efficient Drug Deposition
- ~ Higher Dosing Accuracy
- ~ Portable

## DEVICE CHARACTERISTICS

Operating Humidity Range Max(%rh)	-	95
Operating Humidity Range Min(%rh)	-	10
Operating Temperature Range Max(Degc)	-	33.8 Degree C
Operating Temperature Range Min(Degc)	-	25 Degree C
Storage Humidity Range Max(%rh)	-	95
Storage Humidity Range Min(%rh)	-	10
Storage Temperature Range Max(Degc)	-	3.8 Degree C
Storage Temperature Range Min(Degc)	-	25 Degree C

## DIMENSIONS

Depth (cm)	-	11 cm
Height (cm)	-	22 cm
Weight (Gram)	-	1.4 g
Width (cm)	-	22 cm

## MEDICAL SPECIFICATIONS

Average Nebulization Rate (ML/min)	-	0.4 ML/min
Medication Capacity (ml)	-	3 ml
Mmad (Micro m)	-	2.44 Micro m
Particle Size (Micro m)	-	0.5 to 6 Micro m

## POWER REQUIREMENTS

Power Consumption	-	50 W
Power Required	-	220v -50 Hz
Power Source	-	Electricity

**Hitwin**  
HEALTHCARE

Hitwin Healthcare Pvt. Ltd.  
7/10, Jaganathan Nagar 2nd Main Road,  
Arumbakkam, Chennai, India - 600106.  
Ph: 91 44 4684 5555 (20 Lines), 91 44 4233 3313  
Mob: 91 93441 66778

