

yuwell
You, We Love Life

Hitwin
HEALTHCARE
Respiratory • Sleep Therapy • Cardiology



Catalog

Electromedical products

yuwell

Leading Distributor

Hitwin Healthcare Pvt. Ltd.

7/10, Jaganathan Nagar 2nd Main Road,
Arumbakkam, Chennai, India - 600106.

Ph: 91 44 4684 5555 (20 Lines), 91 44 4233 3313

Mob: 91 93441 66778

www.hitwinhealthcare.com

*Waiting to hear
from you...*



About yuwell

China's largest homecare medical equipment listed company

Yuwell-Jiangsu Yuyue medical equipment & supply Co., Ltd. is a Chinese listed company founded in 1998. Yuwell brings professional health management concept and advanced product solution into daily life, makes a health ecosystem consisting of homecare medical, clinic medical and internet medical and builds up a professional and comprehensive medical service platform.

Yuwell's headquarter is located in Shanghai, China. It owns seven R&D centers and five production bases, which are spread in San Diego(US), Tuttlingen(German), Taiwan, Beijing, Shanghai, Nanjing, Suzhou and Danyang. Besides, its representative offices site over the world, reaching a complete network of R&D, production, sales and service.

Yuwell now is cooperating with Chinese famous medical agencies and core expert group to meet and exceed healthcare requirements in homecare field and chronic disease management field. Seven categories of Yuwell's market share have kept No.1 in three fields of respiratory system, cardiovascular system and endocrine system for many years. Undoubtedly, Yuwell is the top one brand in Chinese homecare medical device industry of 2015.

As internet is developing fast all over the world, Yuwell grasps the high-end technology between internet and medical service. In early 2015, Yuwell set up the Suzhou Cloudoc Health Management Company, aiming at remote chronic diseases management and diagnostic service. With cloud technology and analysis of big data, Yuwell now is building industry 4.0 and smart Yuwell to lead the new era of medical device industry.

In the future, Yuwell will continuously adopt first-class technology in the world and bring it into application. With enthusiasm and hope, all Yuwellers will keep healthcare heart, dedicating into life science research and application. Yuwell. You, we love life.




The Global Marketing Network



Yuyue medical has been devoting to further expansion globally. We are now selling our products in more than 110 countries and regions, and gain good reputation from our business partners.

Certification

- In 2000, wheelchair won American FDA510K certificate.
- In 2002, wheelchair won Japanese SG certificate.
- In 2004, oxygen regulator passed American WHA test.
- In 2005, oxygen concentrator, air-compressing nebulizer and electric suction apparatus won CE certificate.
- In 2006, oxygen concentrator, overbed table, oxygen regulator, ultrasonic nebulizer won American FDA registration.
- In 2008, ISO 13485:2003 quality system certificate.
- In 2008, ISO 9001:2000 quality system certificate.
- In 2009, high frequency medical diagnostic X-ray machine through the China National Compulsory Product Certification CCC Certification.
- In 2009, oxygen concentrator won American FDA510K certificate.
- In 2011, the first won "first prize of excellent quality management team for the pharmaceutical industry" medical device manufacturing enterprises.
- In 2011, oxygen concentrator, electronic blood pressure monitor, wheelchair, electric suction apparatus won "International Standard Product Marking Certificate" .
- In 2011, electronic blood pressure monitor and portable phlegm suction unit won European CE certificate.



Group Brands



Social Responsibilities

Yuyue medical has always been enthusiastic about social responsibilities. In 2014, yuyue medical sponsored oxygen concentrators and air purifiers for the 2nd Youth Olympic Games held in Nanjing. And in the World Expo 2010 Shanghai, Yuyue medical was the medical equipment and service supplier for the Theme Pavilions. There were 2,000 Yuwell wheelchairs donated to China Disabled Person Association.



yuwell Nebulizer



403H Air-compressing Nebulizer

It adopts oil-free high-efficiency valve without the need for daily lubrication and it has low noise level. It can break down liquid into tiny particles similar to gas with extremely strong air compression technology.

It is an ideal product for family and medical professionals, suitable for all kinds of patients.

Power supply: AC220 ± 22V, 50 ± 1Hz
 Input power: ≤ 150VA
 Max negative pressure: ≥ 0.15MPa
 Free air flux compressing pump: ≥ 10L/min
 Max nebulizing rate: ≥ 0.2mL/min
 Operation noise: ≤ 55dB(A)
 Dimension: 19.8cm x 19.5cm x 10.5cm
 Weight: 2.1kgs
 Mass median diameter(MMD): 3.9 μm Percent of MMD < 5 μm: > 65%※
 Work system: Intermittently running & continuous running (24hours)



405A/405B Air-compressing Nebulizer

Input: AC100-220V, 50/60Hz 350mA
 Output: DC12V, 1000mA
 Max nebulizing rate: ≥ 0.2mL/min
 Mass median diameter(MMD): 3.9 μm
 Percent of MMD < 5 μm: > 65%
 Max negative pressure: ≥ 0.1MPa
 Free air flux: ≥ 3.5L/min
 Noise: ≤ 55dB(A)
 N.W.: ≤ 260g (405A) ; ≤ 250g (405B)
 Dimension: 13.5*6.5*6.8cm (405A) ; 14.5*8*5.6cm (405B)
 Work system: Intermittently running.
 (Running for 20 minutes and stop for 40 minutes.)



HL100A Mesh Nebulizer

Features:
 Dimension: 93x41.5x45mm
 Weight: 85g(without battery)
 Operation noise: ≥ 0.2mL/min
 Max nebulizing rate: ≤ 50dB(A)
 Continuous working time: ≥ 30min



NM211C Mesh Nebulizer

Features:
 Dimension: 126x52.5x52mm
 Weight: 110g(without battery)
 Operation noise: ≥ 0.2mL/min
 Max nebulizing rate: ≤ 50dB(A)
 Continuous working time: ≥ 30min

yuwell Nebulizer



403M Air-compressing Nebulizer

This is a new nebulizer model that converts liquid medication into the state of gas to be inhaled by the patients with the newly air compression technology. It's suitable for all kinds of patients and is an ideal product for families and medical professionals.

Power supply: AC220 ± 22V, 50 ± 1Hz
 Max negative pressure compressing pump: ≥0.15MPa
 Free air flux compressing pump: ≥10L/min
 Max nebulizing rate: ≥0.1mL/min
 Noise: ≤65dB(A)
 Patent design for the clamping groove



Powerful heart

Adopts oil free motor
 Long service life
 Small particles



To cure the following disease effectively:

A cold cough, Bronchial asthma, COPD,
 Pneumonia and Upper respiratory tract infection



Features



One-tap operation

Efficient radiator



Safety hook for atomizing tube

Precision connecting port

403T Air-compressing Nebulizer



403E Air-compressing Nebulizer



403K Air-compressing Nebulizer



Technical data:
 Power supply: AC220 ± 22V, 50 ± 1Hz
 Max negative pressure compressing pump: ≥0.15MPa
 Free air flux compressing pump: ≥10L/min
 Max nebulizing rate: ≥0.1mL/min
 Noise: ≤65dB(A)
 Patent design for the clamping groove

Nebulizer inhaler



BW-4 Mouth type



BW-5 Mask type for adult



BW-6 Mask type for child

This product is suitable for the patient with respiratory disease to get the inhaling cure of the liquid or suspension medicine. It should be used with 403 series air-compressing nebulizer or oxygen concentrator.

FEATURES:

Special two-valve design to ensure enough air and to avoid wastage.
 Special air-cushion design to reduce the influence on nebulize from the untreated air, it can also make atomizing particles to be more stable and evenly spread.
 The material of the inhaler is avirulent PP plastic.
 Use it separately to avoid cross infection.
 Special handle design, much more convenient for handle.

yuwell Nebulizer

403AI Air-compressing Nebulizer

It adopts oil-free high-efficiency valve without the need for daily lubrication and it has low noise level. It can break down liquid into tiny particles similar to gas with extremely strong air compression technology. It is an ideal product for family and medical professionals, suitable for all kinds of patients.

Power supply: AC220 ± 22V, 50 ± 1Hz
 Max negative pressure compressing pump: ≥ 0.15MPa
 Free air flux compressing pump: ≥ 10L/min
 Max nebulizing rate: ≥ 0.1mL/min
 Operation noise: ≤ 65dB(A)
 Patent design for the clamping groove



403B Air-compressing Nebulizer




403C Air-compressing Nebulizer

402A Ultrasonic Nebulizer

Power supply: AC220 ± 22V
 Ultrasonic frequency: 1.7MHz ± 10%
 Input power: 60VA
 Maximum nebulizing rate: ≥ 3mL/min
 Timing range: 0-60min no level for adjustment
 Continuous loading time: ≥ 4h
 Max medication cup capacity: 150mL
 Max medication capacity for big cup: 350mL

402B Ultrasonic Nebulizer

402AI Ultrasonic Nebulizer



Functions

