

Oxygen Concentrator





8F-5AW Oxygen Concentrator





Medical grade





5L oxygen flow

Flow Range: 0.5 ~ 5L/min Oxygen Concentration: 95.5%-87% Maximum Outlet Pressure: 40 ~ 70kPa Sound Level: Front side 48dB(A)(average), Whole unit 52dB(A) Input Power: 400VA Weight: 15.5kg (34.1 lbs) Dimension: H50cm(19.7") × W39cm(15.4") × D 24.5cm(9.6")

Minimum Operating Time: 30minutes





Smaller


User-preferred neutral color&size to help blend in to any home environment





Ultra-quiet

Lower than lower noise. Good night with peace and love using yuwell Oxygen Concentrators.



Easily move, light in weight, easy user interface.







Nebulizer therapy function

SPEAK DIRECTLY TO US



MARKETING BY

Hitwin Healthcare Pvt. Ltd.

7/10, Jaganathan Nagar 2nd Main Road, Arumbakkam, Chennai, India - 600106.

Ph: 91 44 4684 5555 (20 Lines), 91 44 4233 3313

Mob: 91 93441 66778

