

KARMA

Karma BT10 Bariatric Wheelchair

**Heavy duty but still lightweight and manoeuvrable,
the BT10 is designed from the ground up to give
mobility to even the most demanding user**



WE
ARE THEIR
WIFE LIFE
MATTERS...
Hitwin
HEALTHCARE

**Self
Propelling
Bariatric
Wheelchair**

Hitwin
HEALTHCARE

Respiratory • Sleep Therapy • Cardiology

Hitwin Healthcare Pvt. Ltd.

7/10, Jaganathan Nagar 2nd Main Road,
Arumbakkam, Chennai, India - 600106.

Ph: 91 44 4684 5555 (20 Lines), 91 44 4233 3313
Mob: 91 93441 66778

www.hitwinhealthcare.com



+ A SIMPLE SOLUTION

Designed for larger clients, the BT10 is an easy to manoeuvre wheelchair that can support users of up to **250 kg**. The adjustable axle position allows for the wheelchair to be set up for the specific needs of the user and allow them to self propel easily. A cost effective prescription wheelchair, that can be built to meet the individual user's needs, but still allow for growth.

Its lightweight aluminium frame is easily folded so it can be lifted and stowed in a vehicle if required. The quick release rear axle allows for the wheels to be removed to further reducing the weight and storage space required.

+ FEATURES

- Angle adjustable backrest (from -5° to +15°)
- Fully detachable and flip back armrests
- Ergonomically designed armrests
- Press down, swingaway detachable legrests with angle adjustable footplates
- Rigidising bar to add strength to frame and to assist the attendant to push the wheelchair if required
- Quick release rear axle
- Swingaway anti-tippers
- Hip belt
- Calf strap
- Dual brake system – attendant brake as well as manual brakes on each wheel
- Pump and tool kit for ongoing standard maintenance



+ ACCESSORIES

- 20" rear wheels
- 200 x 50 front castors
- Height adjustable armrests
- Elevating legrests
- Heel loops
- IV Pole
- Oxygen tank holder

+ SPECIFICATIONS



Seat widths
22", 24", 26"



Seat depths
16", 18", 20"




Rear wheel
self propelling 24"
options: 20"



Front castor
7" x 1 3/4" solid
options:
8" x 2" pneumatic




Seat to floor height
17"
options: 16" - 21"



Anti tippers
swingaway



Back height
16" - 18"
adjustable



Overall width (when open)
seat width +8"



Overall length
125 cm
(with 70° hangers)



Total weight
31.8 kg
(including legrests & wheels)



User weight limit
250 kg
(550 lb)